# Vcheck Feline Tnl

FELINE TROPONIN I

For veterinary use only

## INTENDED USE

The Vcheck Feline Tnl is an *in vitro* diagnostic test kit for the quantitative measurement of cardiac troponin I concentration in feline serum. Troponins are part of the contractile apparatus of cardiac muscle tissue. The troponin complex consists of 3 different proteins (troponin T, C, and I). Of the troponins, cardiac troponin I (Tnl) is most specific for cardiac myocyte injury. Cardiac injury leads to the release of cardiac Tnl into circulation, where its concentration correlates with the degree of cardiac damage. The BIONOTE Vcheck Feline Tnl is designed to be used only by veterinarians.

## PRINCIPLE

The Vcheck Feline Tnl Test Kit is a fluorescent immunoassay for the quantitative measurement of Feline Tnl concentration. The Vcheck Feline Tnl Test Kit uses specific Feline Tnl antibodies that bind to Feline Tnl. The test procedure involves dissolving fluorescent microparticles conjugated to these specific antibodies into assay diluent that reacts with Feline Tnl in the sample. When the sample is applied to the sample hole of the test device, Feline Tnl in the sample migrates along the nitrocellulose membrane and forms complexes with the anti-Feline Tnl antibodies coated on the membrane. As a result, the density of the test line reflects the concentration of Feline Tnl in the sample. The BIONOTE Vcheck Analyzer reads the density of this test line and calculates the Feline Tnl concentration from the calibration curve data. The control line is a reference line that indicates whether the test has been performed correctly.

## MATERIALS PROVIDED

Reagent	5 Tests/Kit
① Vcheck Feline Tnl Test device	5
② Assay diluent tube	5
③ Disposable pipette tip	5
④ Instructions for use	1

## MATERIALS REQUIRED, BUT NOT PROVIDED

- BIONOTE Vcheck Analyzer
- 100 µl pipette

## STORAGE AND STABILITY

- Store the test kit at 1~30 °C. **DO NOT FREEZE.**
- Do not store the test kit in the direct sunlight.
- The test kit is stable until the expiry date that is marked on the package label.

Reagent	Open status	Storage	Stability	Note
Test device	Unopened	1~30 °C, Sealed	12 months	Finished product
	Opened	Do not store	-	Use directly
Assay diluent	Unopened	1~30 °C, Sealed	12 months	Finished product
	Opened	Do not store	-	Use directly

## PRECAUTIONS

- The test kit is for feline use only. Do not use for other animals.
- The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the foil pouch.
- Do not reuse test components.
- Do not touch the membrane in the result window of the test device.
- Do not use the test kit beyond the stated expiry date marked on the label.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not mix components from different lot numbers, the components in this kit have been quality control tested as a standard batch unit.

- All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.
- Hemolyzed serum with precipitate such as fibrin may give erroneous results.
- Strictly follow the test procedure including the amount of sample used, as failure to do so may adversely affect test performance and/or produce invalid results.
- Although the Vcheck Feline Tnl Test Kit offers simple and quick quantitative measurement of Tnl concentration in feline serum, there may be a difference in the detection performance with other clinical or laboratory methods with more sophisticated principles.
- Reference range may vary depending on the test population, and each laboratory shall set its own reference range.
- Professional veterinarian should make a final diagnosis based on the results of this product, other test results and clinical findings.
- BIONOTE Vcheck analyzer is recommended to use at 15~30°C.

## COLLECTION AND PREPARATION OF SAMPLE

- Feline serum should be used with this test.
- [Serum]**
  - Collect the whole blood into a blood collection tube that does not contain anticoagulant. Leave to settle for 30 minutes for blood coagulation and then centrifuge to get a serum supernatant.
  - The samples should be tested immediately after collection. If serum samples are not tested immediately, they can be stored frozen at -20 °C or below for up to 2 months. Samples that have been repeatedly frozen and thawed cannot be used.
  - Serum samples stored frozen at -20 °C or below should be brought to room temperature (15~30 °C) and be completely thawed and mixed via vortexing or pipetting before use.

## TEST PROCEDURE

Allow all kit components and sample to reach room temperature (15~30 °C) prior to testing.

### [Coding]

- Turn on V200 Analyzer and select "Standard Test".


- Remove the test device from the foil pouch. Once the "Insert Device" is displayed in the screen, insert the test device.

### [Dilution of sample & Measurement of Tnl value]


- Use a 100 µl pipette to draw 100 µl of the sample and add to the assay diluent tube (100 µl).

- Use a 100 µl pipette to mix the sample with diluent by pipetting for 5~6 times.


- Add the mixed sample (100 µl) into the sample hole of the test device and press the [START] button to initiate testing.

- \* If the time to press [START] button is delayed, it may affect the test result.

- The V200 analyzer will display the test result on the screen after 10 minutes.

- Remove the test device.

\* **Strictly follow the test procedure including the amount of sample (100 µl) used and the test time (10 min), as failure to do so may adversely affect test performance and/or produce invalid results.**

## INTERPRETATION OF THE RESULT

- Read the concentration value of Feline Tnl appearing on the display of the BIONOTE Vcheck Analyzer. (0.01 ~ 20 ng/mL)
- If "↓ 0.01 ng/mL" appears on the display, it means the concentration of Feline Tnl the specimen is less than 0.01 ng/mL.
- If "↑ 20 ng/mL" appears on the display, it means the concentration of Feline Tnl the specimen is greater than 20 ng/mL.
- If the [Invalid] result appears on the screen, a retest shall be carried out.

## REFERENCE RANGE

< 0.18 ng/mL	0.18~0.28 ng/mL	> 0.28 ng/mL
Normal	Suspected	Abnormal

\* Tnl test results should always be interpreted adjacent to other diagnostic findings such as auscultation, radiography, and echocardiogram.

## SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. <b>Solution:</b> Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. <b>Solution:</b> Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. <b>Solution:</b> Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. <b>Solution:</b> Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. <b>Solution:</b> Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Barcode Error	

Extremely High Total Hemoglobin	The measured total hemoglobin is out of the range of 7 to 23 g/dL. <b>Solution:</b> This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The test is invalid. <b>Solution:</b> Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	The calibration is overdue. <b>Solution:</b> If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Not Supported Device	Loading a test device that is not supported by the analyzer. <b>Solution:</b> Check whether the test device is manufactured by BioNote, Inc.
EEE	Internal error has occurred. <b>Solution:</b> If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.

## PERFORMANCE CHARACTERISTIC

- Measuring Range  
Feline Tnl concentration can be measured within the range of 0.01~20 ng/mL by using a serum. To get numerical results for ↑ 20 ng/mL concentration ranges, samples need to be diluted with negative serum and measured again. To calculate the final Feline Tnl concentration, the measured result must be multiplied by the dilution factor.
- Interference  
No interference was observed for each substance up to the concentration presented in the following table.

Interfering substances	Concentration
Hemoglobin	< 250 mg/dL
Intralipid	< 2,500 mg/dL
Cholesterol	< 200 mg/dL
Bilirubin	< 20 mg/dL
Vitamin C	< 100 mg/dL

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