

Anigen Rapid

# Toxoplasma Ab

## General Information

**Purpose :** Detection feline antibody against *Toxoplasma gondii*

**Principle :** Immunochromatographic assay

**Specimen :** Feline whole blood, plasma or serum

**Materials provided :**

- Anigen Rapid Toxoplasma Ab test device
- Assay diluents bottle
- Anticoagulant tube
- Disposable capillary tube

**Special Features :**

- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity compared with ELISA
  - Sensitivity : 100%
  - Specificity : 98.7%



## Test Procedure

